

Running XDS automatically

USAGE: XDS.py [OPTION]... IMAGE_FILES
FILES is for one or multiple diffraction image files.
OPTIONS:

-h, --help

Print this help message.

-1, -2, -3, -4, -5

Jump directly to a particular step (to avoid rerun initial steps):

- | | |
|------------------------|---------------------------------|
| -1: XYCOOR + INIT | (Image background analysis) |
| -2: COLSPOT | (Collecting spots for indexing) |
| -3: IDXREF | (Auto-indexing) |
| -4: DEFPIX + INTEGRATE | (Integration) |
| -5: CORRECT | (Scaling) |

-a, --anomal

Distinguishes Friedel pairs for scaling, strategy and completeness statistics. Default is no anomalous contribution.

-r, --high-resolution

Set a high resolution cutoff. Default is 0 (no cutoff).

-R, --low-resolution

Set a low resolution cutoff. Default is 9999.

-c, --cell

Set the expected cell.

For example: -c "79 79 38 90 90 90"

-s, --spg

Set the expected space group using either the space group number or simple string.

For example: -s 18 or -s P21212

-i, --xds-input

Give direct XDS Keyword input.

For example: -i "DETECTOR_DISTANCE= 167.0 JOB= IDXREF AIR= 0.002"

-b, --beam-center-optimize

Starting from the initial given values, search and optimize the beam center coordinates (given by -x, -y or extracted from the header).

-S, --strategy

Force to go for calculating strategy (XPLAN) and then stops.

Less used **OPTIONS:**

- Slow,
Set options to process more accurately.
- p, --project
Set the project name. The default is the prefix taken from image names. The working directory will be: xds_process_"project"
- d, --distance
Set the detector to crystal distance.
- O, --oscillation
Set frame oscillation range in degree.
For example: -O 0.5
- x, --beam-x
Set a new value for ORGX: X-coordinates (in pixels) of the detector origin.
- y, --beam-y
Set a new value for ORGY: Y-coordinates (in pixels) of the detector origin.
- v, --verbose
Turn on verbose output.
- w, --wavelength
Set the x-ray wavelength.

USUAL SESSION:

```
> goimg
> XDS.py ../col1_1_*.img
> less xds_process_col1_1/CORRECT.LP (look at the high resolution limit)

> XDS.py -r 1.9 -a -3 -s P3121 -c "59 59 123 90 90 120" ../col1_1_*.img
> cd xds_process_col1_1

> xscale.py XDS_ASCII.HKL ../xds_process_low/XDS_ASCII.HKL
> xdsconv.py XSCALE.HKL 8 Se shelx
> xdsconv.py XSCALE.HKL 8 Se solve
> xdsconv.py XSCALE.HKL 8 Se ccp4 shelx/XDS_ASCII_F4.hkl
> xdsconv.py XSCALE.HKL 8 Se phaser ccp4/XDS_ASCII.mtz
> xdsconv.py XSCALE.HKL cns
```

AUTHOR Pierre Legrand (at) synchrotron-soleil.fr (Proxima1 Beamline)